



Intersurgical AquaFlow™

Bubble humidifier range



Oxygen and Aerosol Therapy ▪ Humidification

AquaFlow™ oxygen bubble humidifiers

A comprehensive range of devices for delivering unheated humidification to your patients. Designed to accept either refillable or Baxter® pre-filled sterile water bottles, the AquaFlow range of bubble humidifiers are suitable for patients receiving variable oxygen therapy.



Code	Description	Box Qty.
1506000	AquaFlow, oxygen bubble humidifier without bottle	30
1507000	AquaFlow, oxygen bubble humidifier with bottle, 500ml	20
1505000	AquaFlow, oxygen bubble humidifier with bottle plus M12 adaptor, 500ml	20
1521000	AquaFlow, homecare oxygen bubble humidifier with bottle, 120ml	20
1520000	AquaFlow, homecare bubble humidifier with bottle plus M12 adaptor, 120ml	20

Make an enquiry

Flowmeter adaptors

Oxygen therapy flowmeters are available with a variety of connection types. The Intersurgical range of adaptors has been designed to permit the attachment of humidification devices with 9/16 UNF Diameter Index Safety System (D.I.S.S) threaded couplings to a variety of flowmeter types.

Please download the explanatory poster entitled Oxygen Therapy Flowmeters to choose the most suitable adaptor for your application:



www.intersurgical.com/products/oxygen-and-aerosol-therapy/flowmeter-adaptors

Make an enquiry

Code	Description	Box Qty.
1508001	Sure Loc™ flowmeter adaptor	45
1506001	BOC® flowmeter adaptor	50
1509001	Flowmeter adaptor M12	50
1504001	Oxygen flowmeter nipple D.I.S.S	50
1505001	Oxygen flowmeter nipple M12	50

Baxter is a registered trademark of Baxter International Inc. BOC is a registered trademark of The BOC Group Ltd.



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
 T: +44 (0)118 9656 300 F: +44 (0)118 9656 356 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 9001:2015, ISO 13485:2003 and ISO 14001:2015

Please think before you print
 Save energy and paper.
 If you must print this information sheet please print it double sided.

IS10.10 • Issue 8 04.18