



Intersurgical EcoLite™ paediatric mask range

The lightest and most comfortable mask available



Oxygen and Aerosol Therapy • Paediatric Masks



Quality, innovation and choice

Interact with us



www.intersurgical.com



Intersurgical EcoLite™ paediatric mask range

Welcome to the Comfort Zone

The unique design and materials used to manufacture the Intersurgical EcoLite mask provides improved patient comfort and a reduced environmental impact.

Comfortable for the Patient

Improved patient comfort has been key to the development of the Intersurgical EcoLite mask. The latest manufacturing technology has enabled us to combine two non-PVC materials in the same product, the polypropylene material forming the body of the Intersurgical EcoLite mask is clear, lightweight and rigid enough to maintain the masks shape, whilst the second softer TPE material is utilised in the seal, which is in contact with the patients face.

▶ Video available at www.intersurgical.com



Incurved nose seal

Improving the fit and eliminating the need for a separate nose clip. Designed to prevent oxygen or aerosolised drugs entering the patient's eyes

Soft face seal

For improved patient comfort

A choice of "under chin" or "on chin" positions

Provides a better fit on a wider range of patient face shapes

Comfortable for the Environment

The Intersurgical EcoLite mask is an important part of our Eco range, designed as part of our ongoing focus on sustainable development. Various groups have questioned the use of PVC in medical products and its impact on the environment. The utilisation of alternative materials has resulted in the elimination of PVC from the Intersurgical EcoLite mask, reducing the environmental impact.



www.intersurgical.com/info/eco



Lighter

50% lighter than a PVC mask, more comfortable for the patient

Exhalation ports

shape and position designed to minimise rebreathed carbon dioxide

Clear mask

permits visualisation of the patient's nose and lips

Make an enquiry

Code	Description	Tube length	Box Qty.
1190015	Intersurgical EcoLite, paediatric, medium concentration oxygen mask		42
1196015	Intersurgical EcoLite, paediatric, medium concentration oxygen mask with tube	1.2m	40
1454015	Cirrus2 nebuliser, paediatric, Intersurgical EcoLite mask kit with tube	2.1m	36
1198015	Intersurgical EcoLite, paediatric, aerosol mask		42
1138015	Silente, Intersurgical EcoLite, paediatric, adjustable venturi mask kit with tube	1.8m	40

Completely PVC free options

Code	Description	Tube length	Box Qty.
1454010	Cirrus2 nebuliser, paediatric, Intersurgical EcoLite mask kit with non-PVC tube	2.1m	36
1196010	Intersurgical EcoLite, paediatric, medium concentration oxygen mask with non-PVC tube	2.1m	40

Lower environmental impact product



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
T: +44 (0)118 965 6300 F: +44 (0)118 965 6356 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 9001:2015, ISO 13485:2016 and ISO 14001:2015

Please think before you print
Save energy and paper.
If you must print this information sheet please print it double sided.

IS10.27 • Issue 6 07.19